


Non Contact Vital Signs Monitoring System

Alec Adamski
Rohan Iyengar
Kedar Manishankar
Sai Sathiesh Rajan


Project Description

- Requesting \$174.34 to build prototype
- Samples and analyze analog signal containing respiration and heart rate
- Display important data on computer UI and mobile app



Planned Workflow

- Pre-existing Doppler radar
- Signal transmission via microcontroller
 - User controls sampling rate of analog input
 - Trigger contact system
- Signal Deconstruction
- Display information on computer UI
 - Data is stored in a database
- Mobile app display



Required Technical Specifications

- 16 Bit resolution for ADC
 - External ADC evaluation board
- 2 Sampling Channels
- 1000 samples/second
- UART MCU-computer connection
- SMA female connector on ADC





Signal Processing


- Variable Sampling Rate
- Blind Signal Separation of the Input Signal
 - Independent Component Analysis
 - RELAX Algorithm
- Signal Separation
 - Band Pass Filters
 - Track the individual signals

Database/Mobile App


Computer:



Mobile:



Constraints



- Real time display
- Delay will be introduced in order to sample and transmit data
- Finite amount of data can be stored in the database for free
- Wired transmission used to save time and avoid design expenditures into wireless protocol

Marketing Analysis

- Vital sign monitoring systems range from \$1350 to \$3895 depending on features and robustness
- Our prototype will not be as expensive but will not have the same features




Demonstration

1. A user will be able to create a new account or add their health data to a precreated existing account.
2. Once the user is logged in, the raw signal data will be displayed on the screen. The user can set sampling rate and frequency range to analyze and then click a button to go to start recording.
3. The user views processed signal data like heartbeat and breathing rate on the computer interface to track the current sensor measurement. The transients will also be displayed on this screen.
4. After a brief delay of 1-2 seconds, can use the phone to access processed data with historical trends on the mobile application.

Schedule

Noncontact Vital Sign Monitoring System – Gantt Chart



Questions?